

# Precise application in complex thin-walled anatomy

*Benefit from an all-in-one solution uniting cutting-edge technology with new flexibility.*

*Inject, cut, coagulate – no need for instrument exchange – for efficient workflow and ease of use.*

**TAILOR YOUR OWN TECHNIQUE.**



# Elevating 3rd space endoscopy

**Thin electrode** for precise applications with limited thermal damage

**Coaxial nozzle** for high-pressure injection

**Flexible probe** for smooth passage even in retroflexion

**Small diameter** for adequate suction through the working channel

**Locking mechanism** for constant electrode length

Marketing claim list 0188716, V005

VIO® 3

APC3

ERBEJET® 2



**THE HYBRIDKNIFE® FLEX UNITES ELECTROSURGERY WITH HYDRODISSECTION**

Powered by VIO® 3

## Inject needle-free



Expand the submucosal space without instrument exchange using high-pressure, needle-free, tissue-selective hydrodissection.

## Cut and Coagulate precisely

Achieve high electrosurgical precision with limited thermal damage using the thin electrode, and apply various coagulation techniques with VIO® modes.

## Flex with stability

Expect smooth passage and stable electrode length, even in difficult angled or retroflexed endoscopic positions.



Tubing length	2300 mm
Outer diameter (distal end)	2.6 mm
Outer diameter (tubing)	2.1 mm
Electrode length	1.5 mm / 2.0 mm
Type	I and T
Retro-flexibility	More than 210°
T-Type 1.5 mm	No. 20150-110
T-Type 2.0 mm	No. 20150-111
I-Type 1.5 mm	No. 20150-112
I-Type 2.0 mm	No. 20150-113

**FIND YOUR PRODUCT**